
$AURA^{\mathsf{TM}}$


Premium Home Hospital Bed











M. Marisek



The Aura™ 48" Extra Wide Bed

Aura™ Premium 39" with Comfort™ Mattress



Aura™ Premium Extra Wide 48" with Dream™ Mattress



Aura™ Premium Espresso Woodgrain



Aura™ Platinum Slate Gray Fabric



Shown with two rails and Dream™ Mattress

The Aura™ Platinum Luxury Hospital Bed

Shown with no rails and Comfort™ Mattress

All SonderCare mattresses are designed to work with your bed system to provide optimal comfort and pressure redistribution





Bed Specifications & Dimensions

- 10" Low Height (17" Top of Mattress)
- 32" High Height (39" Top of Mattress)
- Mattress Size: 79"L x 39"/48"W
- Mattress Height: 7" Foam or 8" Air
- Overall Bed Size: 81"L x 44"/52"W
- · Weight Capacity: 500 lbs.
- Multi-Height Assist Rail
- 4-Wheel Central Brake Pedal



Available Bed Accessories

- Overbed Reading Light
- · Underbed Auto Nightlight
- 8" Length Extension Package
- Additional Set of Rails
- Overhead Trapeze Helper Bar
- Extra Large Overbed Table
- Sheet Sets & Adjustable Pillows







