

AURA™

Premium Home Hospital Bed





“ SonderCare has given my mom her independence back. It’s so much easier for her to get out of bed and it’s really comfortable. She’s so proud of her bed. ”

M. Marisek

Bed Positions



BACKREST



KNEES



CHAIR



ZERO GRAVITY



TRENDELENBURG
& REVERSE



HIGH-LOW



One-Button Chair Position for
Comfortable Eating, Reading
or Watching TV



Maximum High Height for
Better Care and Repositioning



Illuminating Hand Control for High
Visibility Even in Dark Settings




Zero Gravity Position to Help
Alleviate Pressure and Back Pain



Multi-Height Assist Rails
for Security and Assistance
During Transfers



Ultra Low Height for
Easy Transfers and
Fall Prevention



Hidden Casters and One
Central Brake Pedal for
Stability and Easy Movement



Select Your Aura™ Bed Size



Aura™ Premium 39"



Aura™ Premium Extra Wide 48"



Aura™ Companion 78"



Aura™ Companion Bed Headboard
Graphite Gray (top) or Silverstone



Aura™ 39" & 48" Beds are Available in Premium
Espresso Woodgrain or Platinum Slate Gray Fabric




Aura™ Platinum 39" Bed Shown with
two rails and Dream™ Mattress

All functions and positions are identical among all sizes of Aura™ beds.
The two sides of the Aura™ Companion bed can be set up to operate together or individually.


Aura™ Platinum fabric package is an upgrade from the Aura™ Premium.
All functions and positions are identical regardless of style and finish selection.

All SonderCare mattresses are designed to work with your bed system to provide optimal comfort and pressure redistribution

Premium Pressure Redistribution Mattresses:




Dream™ Foam Mattress	Signature™ Hybrid Mattress
Bamboo Quilt Top	Copper Infused Quilt Top
Gel Memory Foam Layer	Orthopedic Foam Top Layer
Pressure Reducing All-Foam Design	Hybrid Foam & Pocket Coil Design
Firm Foundation for Support	
Strong Perimeter for Transfers	
Reversible - Soft & Firm Sides	
Available with Zip-On Fluid-Proof Cover	



Other Specialty Mattresses Available:

Comfort™ Mattress includes Memory Foam without Quilt-Top (Non-Reversible)

Air Mattress includes Alternating Pressure / Low-Air-Loss for wound treatment



Bed Specifications & Dimensions

- 10" Low Height (17" Top of Mattress)
- 32" High Height (39" Top of Mattress)
- Mattress Size: 79"L x 39"/48"/78"W
- Mattress Height: 6" - 9"
- Overall Bed Size: 81"L x 44"/52"/85"W
- Weight Capacity: 500 lbs.
- Multi-Height Assist Rail
- 4-Wheel Central Brake Pedal

Available Bed Accessories

- Additional Set of Rails
- 8" Length Extension Package
- Overhead Trapeze Helper Bar
- Rail Organizer and Rail Pads
- Overbed Reading Light *
- Underbed Auto Nightlight *
- Extra Large Overbed Table
- Sheets, Duvet Sets & Pillows

* Not compatible with the Companion Bed